



**MATRIX HEALTH LINE** &



**USA**

[www.matrixhealthline.com](http://www.matrixhealthline.com)

917-776-5410

[www.cemworldusa.com](http://www.cemworldusa.com)

718-648-1437

**Authorized dealers in USA: Yelena Zinchenko & Michael Pizhadze  
Director, A.M. Kozhemyakin**



**Certificate № ROSS RU.AE88.V00235**



## **ACTIVATOR**


Emitter CEM® ACTIVATOR »№ 1 is designed to work with the device CEM® ACTIVATOR "for treatment of liquids and solutions, placed in standard plastic bottles with a capacity 0,33 - 3l.

Emitter-tube is made in such a way that one could tighten the plastic bottle and inverting it upside down, put the bottle cap down on a flat surface. Thanks to this liquid in the bottle directly in contact with the surface of EHF radiator without risk to fill the radiator fluid processed.



1. Emitter-tube is connected to the device CEM® ACTIVATOR "with a connecting cable with standard USB connectors. Since the transmitter is allowed to use the cable type USB A / miniB 5P with wide base of the tail from the connector miniB not more than 10 mm. Length of connecting cable should be no more than 1,8 m.
2. Do not use cables longer than recommended, as this may lead to a decrease in pulse power of the emitter.

3. Before use, it is recommended that part of the radiator, which will be in contact with the liquid to disinfect a flannel cloth (gauze, cotton wool), slightly dampened with ethyl alcohol. Before disinfecting wipe thoroughly overcome.
4. To prevent liquid from getting into the radiator NEVER wipe radiator with abundantly moistened cloth.
5. To avoid damage to the radiator **NEVER** immerse the transmitter in the water and expose it to shock and excessive mechanical stress.
6. When dealing with biologically active substances and chemicals (fungicides, etc.) should take all necessary precautions associated with the exception of accidental spills into the radiator, and with the exception of accidental contact with skin, mucous membranes, the respiratory tract.
7. Use of biologically active substances and chemicals only in accordance with instructions for their application, as set by the manufacturer. Implement when working with biologically active substances and chemicals, personal protective equipment.
8. **Warning Output radiation of the radiator located on the side opposite to the marking rates of**



**the emitter.**

9. When you work with the device may not look to the emitter of the output radiation with the naked eye (at the location of the emitter at a distance of less than 50 cm from the eye) to avoid possible undesirable physiological effects.
10. **Do not connect the device CEM® ACTIVATOR "and radiators with connection cables for any other equipment (such as a PC).**



*For additional information contact authorize dealers of Cem Tech Products in USA*

[www.cemworldusa.com](http://www.cemworldusa.com)

917-776-5410

[cemtechusa@gmail.com](mailto:cemtechusa@gmail.com)

125 Oceana Drive East, 2-B

[www.matrixhealthline.com](http://www.matrixhealthline.com)

718-648-1437

[matrixofhealth@gmail.com](mailto:matrixofhealth@gmail.com)

Brooklyn, New York 11235